



ITGB1BP2 Antibody

Product Code	CSB-PA011882GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKP3
Immunogen	Human ITGB1BP2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	integrin beta 1 binding protein (melusin) 2;ITGB1BP2;CHORDC3;ITGB1BP;MELUSIN;MGC119214;MSTP015 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ITGB1BP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.