



KCNJ13 Antibody

Product Code	CSB-PA012051GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60928
Immunogen	Human KCNJ13
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	potassium inwardly-rectifying channel, subfamily J, member 13;KCNJ13;KIR1.4;KIR7.1;MGC33328;SVD ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	KCNJ13
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.