



KHDRBS3 Antibody

Product Code	CSB-PA012156GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75525
Immunogen	Human KHDRBS3
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	KH domain containing, RNA binding, signal transduction associated 3;KHDRBS3;Etle;SALP;SLM-2;SLM2;T-STAR;TSTAR;etoile ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	KHDRBS3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.