



KIAA0319L Antibody

Product Code	CSB-PA012182GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IZA0
Immunogen	Human KIAA0319L
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	KIAA0319-like; KIAA0319L; KIAA1837; PKD1-like; PP791; RP4-765A10.3 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	KIAA0319L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.