



KIF4A Antibody

Product Code	CSB-PA012338GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95239
Immunogen	Human KIF4A
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	kinesin family member 4A;KIF4A;FLJ12530;FLJ12655;FLJ14204;FLJ20631;HSA271784;KIF4;KIF4-G1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	KIF4A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.