



KIRREL3 Antibody

Product Code	CSB-PA012372GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IZU9
Immunogen	Human KIRREL3
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	kin of IRRE like 3 (Drosophila);KIRREL3;KIAA1867;KIRRE;MGC126824;MGC126850;MRD4;NEP H2;PRO4502 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	KIRREL3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.