



LARP7 Antibody

Product Code	CSB-PA012756GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q4G0J3
Immunogen	Human LARP7
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	La ribonucleoprotein domain family, member 7;LARP7;DKFZp564K112;HDCMA18P;MGC104360;PIP7S ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	LARP7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.