



LGALS3 Antibody

Product Code	CSB-PA012887GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17931
Immunogen	Human GAL3
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	lectin, galactoside-binding, soluble, 3;LGALS3;CBP35;GAL3;GALBP;GALIG;LGALS2;MAC2 ;Galectin 3
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	LGALS3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.