



LPCAT1 Antibody

Product Code	CSB-PA013052GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NF37
Immunogen	Human LPCAT1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	lysophosphatidylcholine acyltransferase 1;LPCAT1;AYTL2;FLJ12443;FLJ41609;PFAAP3;lpcat ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	LPCAT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.