



LPCAT4 Antibody

Product Code	CSB-PA013055GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q643R3
Immunogen	Human LPCAT4
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	lysophosphatidylcholine acyltransferase 4;LPCAT4;AGPAT7;AYTL3;FLJ10257;LPAAT-eta;LPEAT2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	LPCAT4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.