



LRRC8D Antibody

Product Code	CSB-PA013171GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7L1W4
Immunogen	Human LRRC8D
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	leucine rich repeat containing 8 family, member D;LRRC8D;FLJ10470;FLJ20403;LRRC5 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	LRRC8D
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.