



LSM14A Antibody

Product Code	CSB-PA013203GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8ND56
Immunogen	Human LSM14A
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	LSM14A, SCD6 homolog A (S. cerevisiae);LSM14A;C19orf13;DKFZp434D1335;FAM61A;RAP55 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	LSM14A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.