



LSM4 Antibody

Product Code	CSB-PA013207GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y4Z0
Immunogen	Human LSM4
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	LSM4 homolog, U6 small nuclear RNA associated (S. cerevisiae);LSM4;YER112W ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	LSM4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.