



MAGT1 Antibody

Product Code	CSB-PA013368GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H0U3
Immunogen	Human MAGT1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	magnesium transporter 1;MAGT1;DKFZp564K142;FLJ14726;IAP;MGC64926;OST3B;PRO0756;bA217 H1.1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MAGT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.