



MAPK12 Antibody

Product Code	CSB-PA013451GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53778
Immunogen	Human MAPK12
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	mitogen-activated protein kinase 12;MAPK12;ERK3;ERK6;P38GAMMA;PRKM12;SAPK-3;SAPK3 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MAPK12
Target Details	Activation of members of the mitogen-activated protein kinase family is a major mechanism for transduction of extracellular signals. Stress-activated protein kinases are one subclass of MAP kinases. This protein functions as a signal transducer during differentiation of myoblasts to myotubes.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.