



MDM1 Antibody

Product Code	CSB-PA013625GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TC05
Immunogen	Human MDM1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	Mdm1 nuclear protein homolog (mouse);MDM1;FLJ95264 ;
Product Type	Purified Rabbit Anti Mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MDM1
Target Details	This gene encodes a nuclear protein similar to the mouse double minute 1 protein. The mouse gene is located in double minute (DM) chromatin particles and is amplified in the mouse transformed 3T3 cell line, and the protein is able to bind to p53. In mouse several transcripts have been described for this gene which result from alternative polyadenylation, splicing and exon usage.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.