



MEMO1 Antibody

Product Code	CSB-PA013693GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y316
Immunogen	Human MEMO1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	mediator of cell motility 1;MEMO1;C2orf4;CGI-27;DKFZp434I0135;FLJ25031;MEMO;NS5ATP7 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MEMO1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.