



MFAP5 Antibody

Product Code	CSB-PA013750GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13361
Immunogen	Human MFAP5
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC
Form	Liquid
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MFAP5