



CENPU Antibody

Product Code	CSB-PA014622GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q71F23
Immunogen	Human MLF1IP
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	MLF1 interacting protein;MLF1IP;CENP-50;CENP-U;CENP-U(50);CENP50;CENPU;FLJ23468;KLIP1;PBIP1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CENPU
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.