



MLH1 Antibody

Product Code	CSB-PA014624GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P40692
Immunogen	Human MLH1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	mutL homolog 1, colon cancer, nonpolyposis type 2 (E. coli);MLH1;COCA2;FCC2;HNPCC;HNPCC2;MGC5172;hMLH1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MLH1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.