



MLLT4 Antibody

Product Code	CSB-PA014636GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55196
Immunogen	Human MLLT4
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	0.1M NaHCO ₃ ,0.1M Glycine, 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	myeloid/lymphoid or mixed-lineage leukemia (trithorax homolog, Drosophila); translocated to, 4;MLLT4;AF-6;AF6;AFADIN;FLJ34371;RP3-431P23.3 ;
Product Type	Purified Rabbit Anti Human,Mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	AFDN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.