



MLXIP Antibody

Product Code	CSB-PA014644GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HAP2
Immunogen	Human MLXIP
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	MLX interacting protein;MLXIP;KIAA0867;MIR;MONDOA;bHLHe36 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Transcription
Target Names	MLXIP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.