



# MMP28 Antibody

|                            |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA014675GA01HU                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Uniprot No.</b>         | Q9H239                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Immunogen</b>           | Human MMP28                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Raised In</b>           | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Species Reactivity</b>  | Human,Mouse,Rat                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Tested Applications</b> | ELISA,WB                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Storage Buffer</b>      | PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Purification Method</b> | Antigen Affinity purified                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Isotype</b>             | IgG                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Alias</b>               | matrix metalloproteinase 28;MMP28;EPILYSIN;MM28;MMP-28;MMP25 ;                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Product Type</b>        | Purified Rabbit Anti human PolyClonal Antibody                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Target Names</b>        | MMP28                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Target Details</b>      | Proteins of the matrix metalloproteinase (MMP) family are involved in the breakdown of extracellular matrix for both normal physiological processes, such as embryonic development, reproduction and tissue remodeling, and disease processes, such as asthma and metastasis. This gene encodes a secreted enzyme that degrades casein. Its expression pattern suggests that it plays a role in tissue homeostasis and in wound repair. Transcript variants encoding different isoforms have been described. |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.                                                                                                                                                                                                                                                                                                                                                                                                                                  |