



MOB1A Antibody

Product Code	CSB-PA014696GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H8S9
Immunogen	Human MOBKL1B
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	MOB1, Mps One Binder kinase activator-like 1B (yeast);MOBKL1B;C2orf6;FLJ10788;FLJ11595;MATS1;MOB1;MOBK1B;Mob4B ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MOB1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.