



MPG Antibody

Product Code	CSB-PA014749GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29372
Immunogen	Human MPG
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	N-methylpurine-DNA glycosylase;MPG;AAG;APNG;CRA36.1;MDG;Mid1;PIG11;PIG16;anpg ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MPG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.