



MPP6 Antibody

Product Code	CSB-PA014763GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZW5
Immunogen	Human MPP6
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	membrane protein, palmitoylated 6 (MAGUK p55 subfamily member 6);MPP6;PALS2;VAM-1;VAM1;p55T ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MPP6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.