



MPRIP Antibody

Product Code	CSB-PA014769GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6WCQ1
Immunogen	Human MPRIP
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	myosin phosphatase Rho interacting protein;MPRIP;KIAA0864;M-RIP;RHOIP3;p116Rip ;
Product Type	Purified Rabbit Anti Human,Mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MPRIP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.