



MST1R Antibody

Product Code	CSB-PA015056GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q04912
Immunogen	Human MST1R
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	macrophage stimulating 1 receptor (c-met-related tyrosine kinase);MST1R;CD136;CDw136;PTK8;RON ;
Product Type	Purified Rabbit Anti Homo sapiens PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MST1R
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.