



MTHFD1L Antibody

Product Code	CSB-PA015154GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6UB35
Immunogen	Human MTHFD1L
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	methylenetetrahydrofolate dehydrogenase (NADP+ dependent) 1-like;MTHFD1L;DKFZp586G1517;FLJ21145;FTHFSDC1;MTC1THFS;dJ292B18.2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MTHFD1L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.