



MTMR14 Antibody

Product Code	CSB-PA015178GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NCE2
Immunogen	Human MTMR14
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	myotubularin related protein 14;MTMR14;C3orf29;FLJ22405;FLJ46453;FLJ90311;HJUMPY;hEDTP ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MTMR14
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.