



MUL1 Antibody

Product Code	CSB-PA015235GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q969V5
Immunogen	Human MUL1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	mitochondrial E3 ubiquitin ligase 1;MUL1;C1orf166;FLJ12875;GIDE;MAPL;MULAN;RNF218;RP11-401M16.2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MUL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.