



MUS81 Antibody

Product Code	CSB-PA015240GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96NY9
Immunogen	Human MUS81
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3.-20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	MUS81 endonuclease homolog (S. cerevisiae);MUS81;FLJ21012;FLJ44872 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MUS81
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.