



MYLK2 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA015321GA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9H1R3 |
| Immunogen | Human MYLK2 |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA, WB |
| Storage Buffer | PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| Purification Method | Antigen Affinity purified |
| Isotype | IgG |
| Alias | myosin light chain kinase 2; MYLK2; KMLC; MLCK; MLCK2; skMLCK ; |
| Product Type | Purified Rabbit Anti human PolyClonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | MYLK2 |
| Target Details | This gene encodes a myosin light chain kinase, a calcium/calmodulin dependent enzyme, that is exclusively expressed in adult skeletal muscle. |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |