



MYLPF Antibody

Product Code	CSB-PA015326GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96A32
Immunogen	Human MYLPF
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	myosin light chain, phosphorylatable, fast skeletal muscle;MYLPF;DKFZp779C0757;HUMMLC2B;MGC13450;MRLC2;MYL11 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MYLPF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.