



NAPG Antibody

Product Code	CSB-PA015450GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99747
Immunogen	Human NAPG
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	N-ethylmaleimide-sensitive factor attachment protein, gamma; NAPG; GAMMASNAP ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NAPG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.