



NCAPD3 Antibody

Product Code	CSB-PA015515GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42695
Immunogen	Human NCAPD3
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	non-SMC condensin II complex, subunit D3;NCAPD3;CAP-D3;FLJ42888;KIAA0056;MGC104671;hCAP-D3;hHCP-6;hcp-6 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NCAPD3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.