



NDUFA4L2 Antibody

Product Code	CSB-PA015627GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NRX3
Immunogen	Human NDUFA4L2
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB, IHC, IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 4-like 2; NDUFA4L2; FLJ26118; NUOMS ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NDUFA4L2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.