



NDUFA6 Antibody

Product Code	CSB-PA015631GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56556
Immunogen	Human NDUFA6
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 6, 14kDa;NDUFA6;B14;CI-B14;LYRM6;NADHB14 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NDUFA6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.