



NDUFA9 Antibody

Product Code	CSB-PA015634GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16795
Immunogen	Human NDUFA9
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 9, 39kDa;NDUFA9;MGC111043;NDUFS2L;SDR22E1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NDUFA9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.