



NEIL1 Antibody

Product Code	CSB-PA015693GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96FI4
Immunogen	Human NEIL1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	nei endonuclease VIII-like 1 (E. coli);NEIL1;FLJ22402;FPG1;NEI1;hFPG1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NEIL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.