



NFE2L1 Antibody

Product Code	CSB-PA015751GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14494
Immunogen	Human NFE2L1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	nuclear factor (erythroid-derived 2)-like 1;NFE2L1;FLJ00380;LCR-F1;NRF1;TCF11 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NFE2L1
Target Details	This gene encodes a protein that is involved in globin gene expression in erythrocytes. Confusion has occurred in bibliographic databases due to the shared symbol of NRF1 for this gene, NFE2L1, and for nuclear respiratory factor 1 which has an official symbol of NRF1.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.