



NHP2L1 Antibody

Product Code	CSB-PA015794GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55769
Immunogen	Human NHP2L1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	NHP2 non-histone chromosome protein 2-like 1 (S. cerevisiae);NHP2L1;15.5K;FA-1;FA1;NHPX;OTK27;SNRNP15-5;SNU13;SPAG 12;SSFA1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SNU13
Target Details	Originally named because of its sequence similarity to the Saccharomyces cerevisiae NHP2 (non-histone protein 2), this protein appears to be a highly conserved nuclear protein that is a component of the [U4/U6.U5] tri-snRNP. It binds to the 5 stem-loop of U4 snRNA. Two transcript variants encoding the same protein have been found for this gene.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.