



NIPSNAP3B Antibody

Product Code	CSB-PA015820GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BS92
Immunogen	Human NIPSNAP3B
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	nipsnap homolog 3B (C. elegans); NIPSNAP3B; FLJ11275; NIPSNAP3 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NIPSNAP3B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.