



NLGN4Y Antibody

Product Code	CSB-PA015858GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NFZ3
Immunogen	Human NLGN4Y
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	neuroligin 4, Y-linked; NLGN4Y; KIAA0951 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NLGN4Y
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.