



NOL8 Antibody

Product Code	CSB-PA015926GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q76FK4
Immunogen	Human NOL8
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	nucleolar protein 8;NOL8;C9orf34;DKFZp686P12242;FLJ20736;Nop132;bA62C3.3;bA62C3.4 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NOL8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.