



NUBP2 Antibody

Product Code	CSB-PA016143GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5Y2
Immunogen	Human NUBP2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	nucleotide binding protein 2 (MinD homolog, E. coli);NUBP2;CFD1;NUBP1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NUBP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.