



OIP5 Antibody

Product Code	CSB-PA016316GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43482
Immunogen	Human OIP5
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	Opa interacting protein 5; OIP5; 5730547N13Rik; LINT-25; MIS18beta ; Protein Mis18-beta;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	OIP5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.