



OLFM3 Antibody

Product Code	CSB-PA016322GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96PB7
Immunogen	Human OLFM3
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	olfactomedin 3;OLFM3;NOE3;NOELIN3;NOELIN3_V1;NOELIN3_V2;NOELIN3_V3;NOELIN3_V4;NOELIN3_V5;NOELIN3_V6;OPTIMEDIN ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	OLFM3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.