



# PARD6B Antibody

|                            |                                                                                                                                                                                                                                                                                                                            |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA017447GA01HU                                                                                                                                                                                                                                                                                                         |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                                                                                                                                                                                                                              |
| <b>Uniprot No.</b>         | Q9BYG5                                                                                                                                                                                                                                                                                                                     |
| <b>Immunogen</b>           | Human PARD6B                                                                                                                                                                                                                                                                                                               |
| <b>Raised In</b>           | Rabbit                                                                                                                                                                                                                                                                                                                     |
| <b>Species Reactivity</b>  | Human,Mouse,Rat                                                                                                                                                                                                                                                                                                            |
| <b>Tested Applications</b> | ELISA,WB,IHC,IF                                                                                                                                                                                                                                                                                                            |
| <b>Storage Buffer</b>      | PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.                                                                                                                                                                                                                                      |
| <b>Purification Method</b> | Antigen Affinity purified                                                                                                                                                                                                                                                                                                  |
| <b>Isotype</b>             | IgG                                                                                                                                                                                                                                                                                                                        |
| <b>Alias</b>               | par-6 partitioning defective 6 homolog beta (C. elegans);PARD6B;PAR6B ;                                                                                                                                                                                                                                                    |
| <b>Product Type</b>        | Purified Rabbit Anti human PolyClonal Antibody                                                                                                                                                                                                                                                                             |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                                                                                                                                                                                                                                                                       |
| <b>Target Names</b>        | PARD6B                                                                                                                                                                                                                                                                                                                     |
| <b>Target Details</b>      | This gene is a member of the PAR6 family and encodes a protein with a PSD95/Discs-large/ZO1 (PDZ) domain, an OPR domain and a semi-Cdc42/Rac interactive binding (CRIB) domain. This cytoplasmic protein is involved in asymmetrical cell division and cell polarization processes as a member of a multi-protein complex. |