



PARP15 Antibody

Product Code	CSB-PA017462GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q460N3
Immunogen	Human PARP15
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	poly (ADP-ribose) polymerase family, member 15; PARP15; BAL3; FLJ40196; FLJ40597; MGC126750; MGC126752 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PARP15
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.